

# PeptiQuant™ Plus

Instrument Performance Kit

Quality Control Kits





# Quality Control Kit

Instrument Performance Kit

## Quality Control Kits

### Everyday Quality Control for Instrument Performance Testing

- Instrument Performance Kits offer everyday quality assurance in the form of a ready-to-use standard mix for routine LC-MS instrument monitoring
- Specifically designed for quality control of LC-MS based proteomics workflows
- Kits include peptides spanning a range of masses and chemistries to check performance on a wide variety of analytes
- Fully-optimized for a variety of LC-MS platforms from different vendors



## PEPTIQUANT™ PLUS QUALITY CONTROL KITS


At MRM Proteomics Inc., we understand the importance of accurate and reproducible quantitation. That's why we developed PeptiQuant™ Plus Quality Control Kits, which are specially optimized tools to assess the performance of quantitative MRM-MS proteomics workflows. PeptiQuant™ Plus Quality Control Kits are offered in two types: Instrument Performance Kits and Workflow Performance Kits.

**Instrument Performance Kits** offer everyday quality assurance in the form of a ready-to-use standard mix that can be quickly re-suspended and injected for routine instrument monitoring. This kit makes it easy to recognize retention time drift or changes in instrument

performance. Because they are specifically designed to test instrument performance for proteomics analyses, they can identify issues such as loss of sensitivity that would be missed using common calibration standards.

Each PeptiQuant™ Plus Quality Control Kit comes complete with internal standards, instrument files, quantitation software and clear instructions. Whether you are running your first PeptiQuant™ Plus Proteomics Kit or you have been running your own MRM-MS proteomics assays for decades, our QC Kits are there to help ensure that your proteomics workflow and instrumentation is up to standard for every sample you analyze.

Instrument Performance Kit reveals differences in instrument sensitivity over time and between instruments





## **MRM** Proteomics

MRM Proteomics  
SB-5100  
141 Ave du Président-Kennedy  
Montreal, QC, H2X 1Y4 Canada

+1 (800) 559-9921  
[info@mrmproteomics.com](mailto:info@mrmproteomics.com)  
[mrmproteomics.com](http://mrmproteomics.com)